Product name	ISO13485 Factory Wholesales 3P Metal Lemo Male to 4.0 Electrode Female Snap Physiotherapy Lead Wire for Electrode Tens Machine
application	Stimulation machine, physical therapy instrument
length	10cm, 20cm, 50cm, 100cm or customized
color	Black,Gray, white, black and customized
Conductor material	Tinned copper
Package	Individually packaged

#### **Purpose and Application Scenarios**

- **Hospital Physiotherapy Departments:** Widely used in professional medical settings for patients undergoing post surgery rehabilitation, treating sports injuries, or managing chronic pain conditions.
- **Rehabilitation Centers:** Essential for providing effective physiotherapy services to clients recovering from strokes, spinal cord injuries, and other physical disabilities.
- **Home based Physiotherapy:** Convenient for patients who prefer to receive physiotherapy treatments at home under the guidance of a healthcare professional.

#### Main features:

- **3P Metal Lemo Male Connector:** The <u>3 pin metal Lemo male connector cable</u> is highly durable and corrosion resistant. Its metal construction ensures a secure and stable connection, minimizing the risk of signal loss. The precision engineered pins are designed to fit snugly into compatible sockets, providing a reliable electrical connection.
- **4.0 Electrode Female Snap:** The **4.0 electrode female snap cable** is designed for easy attachment and detachment of electrodes. It has a robust locking mechanism that keeps the electrodes firmly in place during treatment, preventing accidental disconnection. The buckle is made of high quality plastic that is both lightweight and resistant to wear and tear.
- 1 to 2 Splitting Design: The unique one to two splitting design allows for the simultaneous connection of two electrodes to a single source, increasing treatment efficiency. Such as the lemo redel connectors 6pin to magnetic snap electrode lead wire. This design is especially useful for treating larger areas of the body or for performing bilateral treatments.





#### **Certification and Compliance:**

This cable meets strict standards to ensure the highest quality and safety levels in the medical field:

- **ISO 13485:** Certification ensures compliance with international standards for medical device quality management systems.
- ISO 9001: Ensure that products are manufactured under a sound quality management system.
- **FDA Approved:** The product is registered with the U.S. Food and Drug Administration and meets its medical device standards.
- **CE MDR certification:** Complies with the EU Medical Device Regulation (MDR) to ensure safety and effectiveness for medical use.

## Company information. Overview



10+ Years Experienced Professional Manufacturer of Customized Cables.





100+ Employee





OEM/ODM



2000 M² Factory space with 10years of production experience

















## **Certificate**

















### **FAQ**

### Q1: Can I start a sample order before ordering more?

A: Of course, Heymed will provide free samples if we have them in stock, customized samples require installation fees. Thank you for asking us for samples to check our quality.

### Q2: What payment methods do you accept?

A: PayPal, T/T, bank transfer, Western Union and L/C

#### Q3: How does your factory do in terms of quality control?

Answer: Quality is priority. The QC department with professional quality assurance experts attaches great importance to quality control from beginning to end. Our factory has obtained Rohs, SGS, CE FCC, ISO9001:2008 certification. ETC

## Q4: Factory? Yes! Can I customize the items I like, such as printing our own logo or modifying the design?

A: Yes, we are a factory with more than 100 employees. We provide OEM/ODM services. Please contact us or our sales representative for details.

# Q5: Delivery time and warranty? What shipping service will you use to ship my product?

A: Our delivery time is about 2-4 weeks, and we can provide 12 months warranty from delivery. We use DHL, Fedex, UPS, TNT, EMS as customers.

#### Q6: Where is your factory located? How can I visit there?

Answer: Of course. We are very happy to welcome you to Dongguan and show you our workplace. Our factory is located in the third. Building A, No. 8, Kangyang 1st Street, Shebei Village, Huangjiang Town, Dongguan City, Guangdong Province, China