

品名	4Pin Lemo Lemo 316 ISO13485 DongGuan
规格	
数量	100002000050000100000000
单位	
品牌	
产地	
备注	

4pin Lemo Lemo ISO13485 316 EMS 4PIN Lemo PH2.0-4PIN ISO13485

Lemo 316 Lemo 6PIN DB9SPO2 Lemo

TPU OEM ODM





特長

医療機器に適合するケーブルです。

- **ISO 13485** 認証取得済み
- **ISO 9001** 認証取得済み
- **FDA** 承認済み
- **CE MDR** 認証取得済み

仕様



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



公司

Q1 公司
Heymed

Q2
PayPal T/T L/C

Q3
QC ROHS SGS CE FCC ISO9001 2008 ETC

Q4
100 OEM/ODM

Q5
2-4 12 DHL FedEx UPS TNT EMS

Q6
Dongguan 8 Kangyang 1st Street