

	4pin lemo 40mm 4 Dongguan
	Lemo 4 ECG
	1M
	EEG ECG EMG
	100pcs

ECG/EKG Lemo ECG EEG EEG EKG EMG
 Hollter ECG Snap Cable TPU

EEG EEG EGG EMG





□□□□□□



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



Q1

Heymed

Q2□□□□□□□□□□

PayPal T/T L/C

Q3

QC ROHS SGS CE FCC ISO9001 2008 ETC

4

100 OEM/ODM

5

2-4 12 DHL FedEx UPS TNT EMS

6

[Dongguan](#)[8](#)[1](#)[Shebei](#)
[Village](#)[Shebei](#)