

產品

產品	2PIN LEMO 2.0mm Dongguan ISO13485
規格	2.0mm
長度	1M
顏色	黑色
用途	EEG/ECG/EMG
包裝	100pcs

2PIN LEMO 2.0mm Dongguan ISO13485 TENS/EMS/ECG/EMG

Lemo 2.0mm TPU/PVC

DC2.35mm 3.5mm TENS/EMS 30

USB -C 4.0 earhook ECG/EMG

Dongguan ISO13485

產品





□□□□□□□□



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



ISO13485

FCC

CE

ROHS



UKCA

EPR

EPR

ISO9001

□ □ □ □

Q1

Heymed

Q2

PayPal T/T L/C

Q3□□□□□□□□□□□□□□□□

QC ROHS SGS CE FCC ISO9001 2008 ETC

4

100 OEM/ODM

5

2-4 12 DHL FedEx UPS TNT EMS

6