


品名	2PIN 2.0mm 2.0mm Tens EMS
規格	
長さ	10cm 20cm 50cm 100cm
材質	
色	
単位	
備考	

2PIN 2.0mm 2.0mm Tens EMS

2PIN TENS EMS TPU PVC

DC 2.35mm 2.0mm 2.54mm 2.0mm TENS EMS ECG

OEM/ODM





About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10years of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



ISO 13485

FCC

CE

ROHS



UKCA

EPR

EPR

ISO 9001

