

	5pin Lemo 4.0mm ECG
	Lemo 4.0 ECG
	1M
	EEG ECG EMG
	100pcs


 5PIN LEMO 4.0mm SNAP ECG ECG EMG Lemo
 4.0mm SNAP [USB C 4.0mm](#)
[2](#) [6pin lemo](#)[3.5](#)[2](#)
 ISO13485

 ISO13485 OEM/ODM





□□□□□□□□



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



□ □ □ □

Q1

Heymed

Q2

PayPal T/T L/C

Q3

QC ROHS SGS CE FCC ISO9001 2008 ETC

[illegible]

100 OEM/ODM

5

2-4 12 DHL FedEx UPS TNT EMS

6