

	90C 3 4.0 Dongguan Tens EMS
	3.5mm/4.0mm USB
	120
	TENS EMS / ECG/EKG/EEG/EMG
	500pcs

 90C TO 3 Leads 4.0mm ISO13485 TENS EMS


90 4.0mm TPU PVC

 PIN CES ear

 OEM/ODM n







About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



ISO13485

FCC

CE

ROHS



UKCA

EPR

EPR

ISO9001

□ □ □ □

Q1

Heymed

Q2

PayPal T/T L/C

3

□□□□□□□□□□□□□□QC□□□□□□□□□□□□□□□□□□ROHS□SGS□CE FCC□ISO9001□2008□□□ ETC

4

100 OEM/ODM

5

2-4 12 DHL FedEx UPS TNT EMS

6

□□□□□□□□□□ Dongguan □□□□□□□□□□□□□□□□□□□□ A Building A□8□□□□□□□□□□□□□□ Shebei

□Kangyang 1st Street□Shebei□