

	DB15PIN SR 4.0mm 13
	3.5mm/4.0mm DB 15PIN
	120
	TENS EMS/ ECG/EKG/EEG/EMG
	500pcs

 DB15PIN SR 4.0mm 13

 DB15 SR 10POGO PIN 10
4.0mm ECG

 ECG EMG TENS

 PVC

 2.54 4.0mm

 Dongguan OEM/ODM





□□□□□□



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



ISO13485

FCC

CE

RoHS



UKCA

EPR

EPR

ISO9001

□ □ □ □

Q1

Heymed

Q2□□□□□□□□□□

PayPal T/T L/C

3

□□□□□□□□□□□□□□QC□□□□□□□□□□□□□□□□ROHS□SGS□CE FCC□ISO9001□2008□□□ ETC

[illegible]

100 OEM/ODM

5

2-4 12 DHL FedEx UPS TNT EMS

6

□□□□□□□□□□Dongguan□□□□□□□□□□□□□□□□ A Building A□8□□□□□□□□□□Shebei□

□Kangyang 1st Street□Shebei□