

	Lemo 4 4.0mm Dongguan Tens EMS
	3.5mm/4.0mm Lemo
	120
	TENS EMS / ECG/EKG/EEG/EMG
	500pcs

 Lemo ISO13485 Lemo 4
 TENS EMS ECG EMG

 lemo 4.0mm
TPU PVC


 lemo **4PIN 2.0mm 3**

 OEM/ODM **USB C 4.0mm 2** TENS
EMS C C





□□□□□□



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



□ □ □ □

Q1

Heymed

Q2

PayPal T/T L/C

3

□□□□□□□□□□□□□□QC□□□□□□□□□□□□□□□□□□□ROHS□SGS□CE FCC□ISO9001□2008□□□ ETC

4

100 OEM/ODM

5

2-4 12 DHL FedEx UPS TNT EMS

6

□□□□□□□□□□ Dongguan □□□□□□□□□□□□□□□□□□□□ A Building A□8□□□□□□□□□□□□□□ Shebei

□Kangyang 1st Street□Shebei□