

ISO13485	HEYMED 2.0 TEN 2P PIN -
	10 20 50 100

ISO13485

2.0mm PVC 2.0mm

ISO13485

- 
- PVC
- 1
- 
- 5A

ISO13485

- 
- 
- ECG EMG EEG
- 

ISO13485

ISO13485

ISO13485







## About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M<sup>2</sup>  
Factory space with 10years of production experience



100+  
Employees



OEM/ODM  
Support cable customized



TOP 500  
Bands approved manufacturer



## Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



ISO 13485

FCC

CE

ROHS



UKCA

EPR

EPR

ISO9001

