

	USB 2.0 A 2 Y
	USB A Type C & Micro USB 5pin
	30cm
	100PCS

90YUSB 2.0 A YUSB 2.0 A type c micro usb 5pin


2 USB Y USB Type C .







□□□□□□



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



ISO13485

FCC

CE

ROHS



UKCA

EPR

EPR

ISO9001

--	--	--	--

Q1:

A: Heymed

Q2:

A: PayPal T/T

Q3:

A: 通过 QC 质量管理体系 Rohs SGS CE FCC ISO9001:2008 等 ETC

Q4:

A: 100 OEM/ODM

Q5:

A: 2-4 12 DHL Fedex UPS TNT EMS

Q6:

8A