


	3 4.0 ECG 3.5
	3.5mm4.0mm
	30cm
	100PCS


3 aduio Tens

Heymed ecg eeg ekg emg OEM ODM .
 Micro usb **EEG ECG** .





□□□□□□



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



□□□□

Q1:

A: Heymed

Q2:

A: PayPalT/T

Q3:

[illegible]

Q4:

A: 100 OEM/ODM

Q5:

A: 2-4 12 DHL Fedex UPS TNT EMS

Q6:

8A