

FOB

1,000,000

T/T

3-7

UL

Fedex

500PCS

-
- PSP10
- ED
-

USB C 4.0mm, lemo 2.0mm ETC



00000

000000000000000000000000 3.5mm0000000000000000000000Heymed000000000000
000000000000000000000000

000000 12 00000000 RoHS 0 Reach 00000000

000000I S O 13485000000000000I S O 9 0 0 1000000000I S O 1 4 0 0 1000000000000000
0UL0CQC0CE0SGS000000



11111

Q1:

☐: Heymed ☐

Q2:

□: PayPal□T/T□□□□□□□□□□□□

Q3:

□: □□□□□□□□□□QC□□□□□□□□□□□□□□□Rohs□SGS□CE FCC□ISO9001:2008□□ ETC

Q4:

A: 100 OEM/ODM

Q5:

□: □□□□□□□□ 2-4 □□□□□□□□□□ 12 □□□□□□□□□□ DHL□Fedex□UPS□TNT□EMS □□□□□□

Q6:

[illegible]