


HEYMED







About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



□ □ □ □

Q1

Heymed

Q2

PayPal T/T L/C

Q3

QC ROHS SGS CE FCC ISO9001 2008 ETC

4

100 OEM/ODM

5

2-4 12 DHL FedEx UPS TNT EMS

6