

□□□□

□□□□	8 □ 1 □□□□ □□□□ □□□ 4.0mm □□□□ □□□□ □□
□□□	3.5 □□/4.0 □□□□□□ USB C □□□
□□□□	120cm□□□
□□□□	□□ □□ □□□
□□	TENS EMS □□/□□□ ECG/EKG/EEG/EMG □□□□□
□□□	100PCS

□□□□□□

□□□□□ **USB 3.1 C □□ ECG □□□□**

□□□ 1□USB c □□□□□□□□□□

□□□ 2□4.0mm □□□□□

□□□□□□□□□□□□

□□□□□□□□□□□□


□□ **□□□□C□□□□□□□□□□**□□□□

□□□□□□





□□□□□□



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



□ □ □ □

Q1:

Heymed

Q2:

☐: PayPal

Q3:

□: □□□□□□□□□□QC□□□□□□□□□□□□□□□Rohs□SGS□CE FCC□ISO9001:2008□□ ETC

Q4:

A: 100 OEM/ODM

Q5:

□: □□□□□□□□ 2-4 □□□□□□□□□□ 12 □□□□□□□□□□ DHL□Fedex□UPS□TNT□EMS □□□□

Q6:

8A