



0000	00 LEMO 000000000000
000	00 LEMO 00002.0mm 00
0000	30cm000
0000	00000
0000	000000000000
00	0000
00000	100PCS

4Pin Lemo 2.0mm Din ECG EEG EKG EMG TENS
 Lemo 4 2.0mm Din TPU

□□□□□□□□□□ din □□□□ □□ C □□□□ 2.0mm □□ 4 □□□□ □ .

_____.







About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.





Q1:

A: Heymed

Q2:

A: PayPal T/T

Q3:

A: QC Rohs SGS CE FCC ISO9001:2008 ETC

Q4:

A: 100 OEM/ODM

Q5:

A: 2-4 12 DHL Fedex UPS TNT EMS

Q6:

8A