

	USB 3.0 C 3.5mm
	3.5 /4.0 USB C
	30cm
	TENS EMS ECG/EKG/EEG/EMG
	100PCS

ISO13485

USB C ECG USB C 3.5mm ecg EEG ECG
TPU OEM ODM M5 ECG







About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



--	--	--	--

Q1:

A: Heymed

Q2:

A: PayPal T/T

Q3:

A: QC Rohs SGS CE FCC ISO9001:2008 ETC

Q4:

A: 100 OEM/ODM

Q5:

A: 2-4 12 DHL Fedex UPS TNT EMS

Q6:

8A