

--	--	--	--	--	--	--

厚度	2.54mm
规格	2.54mm 3.5mm/4.0mm ecg
长度	30cm
重量	
材质	TENS EMS ECG/EKG/EEG/EMG
MOQ	100

ISO13485                                                                      <

1. 4.0mm ecg 2.54mm 1x1p 2. OEM ODM 3. RJ45 4. ecg 5. ECG 6. M5

□□□□:







About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



