

Trade Term: FOB Shanghai/Beijing

Quantity: 1,000,000 Pcs

Payment: T/T

Lead Time: 4 weeks

Delivery Address: 370000 15200000

Certification: UL, CE, ROHS, FCC, UKCA


Logistics: Fedex, DHL, UPS

MOQ: 500 PCS

Notes:

1. **Product Description:** The product is a high-precision medical device used for monitoring and recording physiological signals.
2. **Material:** The main body is made of medical-grade stainless steel, ensuring biocompatibility and durability.
3. **Performance:** EEG (Electroencephalogram), ECG (Electrocardiogram), EMG (Electromyography), and TENS (Transcutaneous Electrical Stimulation) signals can be recorded simultaneously with high accuracy.
4. **Lead Time:** 2.0mm lead time: 2.0 mm lead time for standard configurations. Customized configurations may require additional lead time.
5. **3.5mm Lead Time:** 3.5mm lead time for configurations requiring longer leads. Customized configurations may require additional lead time.
6. **Accuracy:** The device is calibrated to ensure high accuracy in signal recording and analysis.
7. **Software:** The device is compatible with our proprietary software for data collection and analysis.
8. **Documentation:** The device is accompanied by detailed user manuals and technical specifications.
9. **Warranty:** The device is covered by a 1-year warranty against manufacturing defects.
10. **After-Sales Service:** We provide comprehensive after-sales support, including technical assistance and spare parts.

3.5mm 连接器 规格书



Heymed Technology 3.5mm 连接器 规格书

12 RoHS Reach

ISO 13485 ISO 9001 ISO 14001 UL CQC CE SGS