

	30cm
	EMG
MOQ	100

 ecg eeg ekg emg ECG DC 2.35
 OEM ODM EEG
ECG .
 USB EEG ECG .

A: 2.2mm

B: 3.5/4.0/5.0mm

 :





□□□□□□



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



ISO13485

FCC

CE

RoHS



UKCA

EPR

EPR

ISO9001

