

	4.0mm 4.0mmPVC TENS EMS Dongguan Tens EMS
	4.0mm
	120cm
	TENS EMS / ECG/EKG/EEG/EMG
	500pcs

 ISO13485 TENS EMS PVC 4.0mm 4.0mm EMS

 2.0mm4.0mmEMS PVC TPU Pin-to -SNAP

 PVC OEM/ODM

 LEMO44.0mmLEMO44.0mm TENS EMS ECG EMG





□□□□□□□□



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



ISO13485

FCC

CE

ROHS



UKCA

EPR

EPR

ISO9001

--	--	--	--	--	--

Q1

A Heymed

Q2

A PayPal/T/T Western Union L/C

Q3

QCRHSGSCE FCC
ISO90012008

Q4

A100 OEM/ODM

Q5

A 24 12 DHL FedEx UPS TNT EMS

Q6

guan3

□A□No□8□Kangyang 1st Street□Shebei Village□Huangjiang Town□Dongguan City□Grina□China