

	4Lemo 40mm 40mm 4
	LEMO 40mm ECG
	1M
	EEG ECG EMG
Moq	100pcs

ECG/EKG Lemo Circular Connector ECG EEG EKG EMG TPU ECG


ECG EEG EKG EMG

.





□□□□□□



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



--	--	--	--	--	--

Q1

A

Q2

A PayPal T/T Western Union L/C

Q3

AQCC QC ROH SGS CE
FCC ISO9001 2008

Q4

A100 OEM/ODM

Q5

A 2412 DHL FedEx UPS TNT EMS

Q6

A 3 A No.8
Kangyang 1st Street Shebei Village Huangjiang Town Dongguan City Grina China