


Moq	500pcs

 Dongguan ISO13485 6 din

 6.35 mm 2 3.5 mm - TENS/EMS Dongguan ISO13485


Precision-Enyered Male Plug PVC

 OEM ODM 316 4 L LEMO TENS EMS ISO 13485





□□□□□□



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



--	--	--	--	--	--

Q1

A

Q2

A PayPal T/T Western Union L/C

Q3

A#####QC#####ROH#SGS#CE
FCC#ISO9001#2008#####

Q4

A100 OEM/ODM

Q5

A 24 12 DHL FedEx UPS TNT EMS

Q6

A No.8
Kangyang 1st Street Shebei Village Huangqiang Town Dongguan City Grina China