

本公司承接各种EEG脑电波信号采集系统BCI脑机接口系统开发及系统集成服务。我们拥有专业的研发团队和丰富的项目经验，能够为客户提供定制化的解决方案。我们的EEG系统广泛应用于医疗、教育、科研等领域，帮助客户实现脑电波信号的精准采集和分析。

我们符合ISO13485和ISO9001质量管理体系标准，确保产品质量。我们提供B2B和B2C服务模式，支持OEM/ODM定制。我们的EEG系统具有高精度、高稳定性、易操作等特点，能够满足不同客户的需求。我们致力于为客户提供优质的产品和服务，提升客户满意度。

我们EEG系统支持多种脑电波信号采集，包括Alpha波、Beta波、Gamma波、Delta波和Theta波。我们提供完整的EEG系统解决方案，包括硬件设备、软件平台和数据分析工具。我们拥有专业的技术支持团队，能够为客户提供及时的技术支持和售后服务。






ISO 13485

ISO 9001

- **ISO 13485** 國際標準化組織 (International Organization for Standardization) 制定之醫療器械品質管理系統標準
- **ISO 9001** 國際標準化組織 (International Organization for Standardization) 制定之品質管理系統標準
- **FDA** 美國食品及藥物管理局 (U.S. Food and Drug Administration) 核准之醫療器械品質管理系統標準
- **CE MDR** 歐盟 (European Union) 核准之醫療器械品質管理系統標準

品質管理系統



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



ISO13485

FCC

CE

RoHS



UKCA

EPR

EPR

ISO9001

Q1
A

Q2
A PayPal T/T Western Union L/C

Q3
 QC ROH SGS CE
FCC ISO9001 2008

Q4
A 100 OEM/ODM

Q5
A 2 4 12 DHL FedEx UPS TNT EMS

Q6
 guan 3
 A No 8 Kangyang 1st Street Shebei Village Huangjiang Town Dongguan City Grina China