

品名	DIN 5 3.5mm ECG guan
規格	
尺寸	10cm 20cm 50cm 100cm
材質	
包裝	
備註	

說明

DIN 5 3.5mm ECG EMG TENS EMS **Spin lemo 4.0mm ECG** **C 4.0mm 2**

DIN 5 2 3.5mm

TPU PVC guan OEM/ODM ISO13485 ROHS







## About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M<sup>2</sup>  
Factory space with 10years of production experience



100+  
Employees



OEM/ODM  
Support cable customized



TOP 500  
Bands approved manufacturer



## Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



ISO 13485

FCC

CE

ROHS



UKCA

EPR

EPR

ISO 9001

