

Table 1

Length	10m/100m/1000m/10000m/4.0mm ECG Cable
Pin Pitch	3.5mm/4.0mm
Cable Length	80cm
Material	PVC
Application	TENS EMS/ECG/EKG/EEG/EMG
Quantity	500pcs

ECG Cable Dongguan ISO 13485 Pogo PIN

4.0mm ECG Cable PVC EEG ECG

9pin 5 2.0mm

Table 2

- 4pin 4.0mm 3 ECG Cable
- 1 4pin Pogo Pin Male Connector
- 2 4.0mm ECG Cable
-
-
- ISO 13485 ISO 9001 CUSOM ECG negative Dongguan ECG Medical Cable Supplier to

Table 3





□□□□□□□□



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



Q1
AHeymed

Q2
APayPalT/TWestern UnionL/C

Q3
 QCROHSGSCE FCC
ISO90012008

Q4
A100 OEM/ODM

Q5
A2412 DHLFedExUPS TNT EMS

Q6
 guan3
ANo8Kangyang 1st StreetShebei VillageHuangjiang TownDongguan CityGrinaChina