

1. 產品規格

產品名稱	2Pin LEMO 2.0mm 2Pin 2Pin Dongguan ISO13485 認證
規格	2Pin LEMO 2.0mm 2Pin
數量	1M
用途	醫療設備
測試項目	EEG ECG EMG
單位	個
Moq	100pcs

1. Dongguan ISO13485 認證，符合醫療設備標準，適用於 2Pin Lemo 2.0mm 2Pin 連接器，適用於 TENS EMS ECG EMG 等醫療設備。

2. LEMO 2.0mm 2Pin 連接器，採用 TPU 材質，具有優良的耐腐蝕性和耐磨性，適用於 PVC 等醫療設備。

3. 產品規格：DC2.35mm 3.5mm TENS/EMS 連接器，長度 30cm。

4. USB - C 連接器，長度 4.0mm，適用於 ECG/EMG 等醫療設備。

5. 產品規格：Dongguan Factory ISO13485 認證，適用於醫療設備。

.

2. 產品規格





□□□□□□□□



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



--	--	--	--	--	--

Q1

A

Q2

A PayPal T/T Western Union L/C

Q3

A#####QC#####ROH#SGS#CE
FCC#ISO9001#2008#####

Q4

A100 OEM/ODM

Q5

A 24 12 DHL FedEx UPS TNT EMS

Q6

A No.8
Kangyang 1st Street Shebei Village Huangqiang Town Dongguan City Grina China