

產品規格:

產品名稱	DC5521 2 in 1 智慧健康監測器 - DC5521 智慧健康監測器
型號	DC5521
尺寸	30cm
重量	1.2kg
功能	TENS EMS 電療/電擊 ECG/EKG/EEG/EMG 心跳監測
電池壽命	500h

特點:

- 產品: DC5521 2 in 1 智慧健康監測器 - DC5521 智慧健康監測器
- 12個月保用: 產品自購買之日起, 保用12個月。如產品在保用期內出現任何故障, 我們將免費為您更換新產品。
- 智慧健康監測: 產品可監測您的心跳、血壓、血氧飽和度等數據, 並可透過智慧手機應用程式查看您的健康狀況。
- 心跳監測: 產品可透過智慧手機應用程式監測您的心跳, 並可透過智慧手機應用程式查看您的心跳數據。
- 電療/電擊: 產品可透過智慧手機應用程式進行電療/電擊, 並可透過智慧手機應用程式查看您的電療/電擊數據。
- 電池壽命: 產品可透過智慧手機應用程式查看您的電池壽命, 並可透過智慧手機應用程式查看您的電池壽命數據。
- 產品: DC5521 2 in 1 智慧健康監測器 - DC5521 智慧健康監測器


產品規格:

- 產品名稱: DC5521 2 in 1 智慧健康監測器 - DC5521 智慧健康監測器
- 型號: DC5521
- 尺寸: 30cm
- 重量: 1.2kg
- 功能: TENS EMS 電療/電擊 ECG/EKG/EEG/EMG 心跳監測
- 電池壽命: 500h

產品規格:

- CE: 符合 ECG 標準
- FDA: FDA 批准 (FDA 批准) 符合 FDA 標準
- RoHS: 符合 RoHS (RoHS) 標準
- PBB: 符合 PBB 標準
- PBDE: 符合 PBDE 標準


產品規格:







□□□□□□



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



ISO13485

FCC

CE

RoHS



UKCA

EPR

EPR

ISO9001

--	--	--	--	--	--

Q1:

A: Heymed

Q2:

A: PayPal T/T L/C

Q3:

QC Rohs SGS CE
FCC ISO9001:2008

Q4

A: 100 OEM/ODM

Q5: □□□□□?□□□□□□□□□□□□□□□□□□□□□?

A: 2 4 12 DHL Fedex UPS TNT EMS

Q6: □□□□□□□□□?□□□□□□□□□□□□□□□□□?

: 38
A