

--	--	--	--	--	--	--

型番	000000000 4PIN Lemo 0000000000 PH2.0-4PIN 0000000000
接続形式	4000000LEMO000000 PH2.0-4PIN00
ケーブル長	2.5M01M01.8M 00000000000000
ケーブル色	000000000000000
測定項目	EEG0ECG0EMG0000
単位	00000000000
MOQ	10

[illegible]

FGG 1B 8Pin LMO - 6Pin
TENS4LEMO

□□□□●●





□□□□□□



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



--	--	--	--	--	--

Q1:

A: Heymed

Q2:

A: PayPal T/T L/C

Q3:

A: ☐ QC ☐ RoHS ☐ SGS ☐ CE
☐ FCC ☐ ISO9001:2008

Q4

A: 100 OEM/ODM

Q5: □□□□□?□□□□□□□□□□□□□□□□□□□□?

A: 2 4 12 DHL Fedex UPS TNT EMS

Q6:

A: 38