

1. 产品概述

规格	2.54mm间距4.0mm间距 Dongguan ISO13485
尺寸	3.5mm/4.0mm间距2.54mm间距
长度	50cm
重量	
应用	TENS EMS/ECG/EKG/EEG/EMG
Moq	500pcs

Dupont 2.54 mm间距4.0 mm间距SNAP PVC 5.0 mm ECG

符合ISO13485 4.0mm

2. 规格





□□□□□□□□



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



ISO13485

FCC

CE

RoHS



UKCA

EPR

EPR

ISO9001

--	--	--	--	--	--

Q1

A

Q2

A PayPal T/T Western Union L/C

Q3

A#####QC#####ROH#SGS#CE
FCC#ISO9001#2008#####

Q4

A100 OEM/ODM

Q5

A 24 12 DHL FedEx UPS TNT EMS

Q6

A No.8
Kangyang 1st Street Shebei Village Huangqiang Town Dongguan City Grina China