

品牌	LEMO 4.0mm Dongguan Medical Wire
规格	
长度	10 cm 20 cm 50 cm 100 cm
颜色	
用途	
其他	

Dongguan ISO 13485 LEMO 4.0mm SANP 4 4.0mm LEMO

4 2.0mm 3

4.0mm 4.0mm TENS EMS





品質保証

品質保証

- **ISO 13485** 品質保証システム認証取得済み
- **ISO 9001** 品質保証システム認証取得済み
- **FDA** 米国FDA承認済み
- **CE MDR** 欧州CE MDR承認済み

品質保証



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



