

	4.0mmRJ45
	RJ45 8P8C 4.0mm
	30cm
	EEGECGEMGTENS
MOQ	100

 RJ45 8P8C 4.0mm ecg EMS TEN

Heymed OEM ODM EEG ECG. 10mmEU.

 :





□□□□□□



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



