

品名	USB Type-C 4mm Crocodile Clips Dongguan
規格	
尺寸	10cm 20cm 50cm 100cm
顏色	
包裝	
材料	
產地	

USB Type-C 4mm Crocodile Clips Dongguan

OEM/ODM USB C 4.0mm - Type - C Dongguan PVC 0.5m ISO13485

4 - in - 1 audiojack 4 ECG 1M TPU 3.5mm 4.0mm Dongguan Factory ISO ISO13485


USB Type - C

Dongguan Factory

-
-
-
-



--	--	--	--	--	--	--	--



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



Q1
A

Q2
A PayPal T/T Western Union L/C

Q3
 QC ROH SGS CE
FCC ISO9001 2008

Q4
A 100 OEM/ODM

Q5
A 2 4 12 DHL FedEx UPS TNT EMS

Q6
 guan 3
 A No 8 Kangyang 1st Street Shebei Village Huangjiang Town Dongguan City Grina China