

Ag/AgCl Electrode Array for ECG, EEG, and EMG Applications

High-quality, disposable Ag/AgCl electrodes for medical monitoring. Available in various sizes and configurations. [View Details](#) [Contact Us](#)

Manufactured by Dongguan Factory, B2B OEM/ODM services available.



□ □ □ □ □ □ □ □



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



--	--	--	--	--	--

Q1

A

Q2□□□□□□□□□□□□□□□□

A PayPal/T/T Western Union L/C

Q3

QC ROH SGS CE
FCC ISO9001 2008

Q4

A100 OEM/ODM

Q5

A 2412 DHL FedEx UPS TNT EMS

Q6

guan3

□A□No□8□Kangyang 1st Street□Shebei Village□Huangjiang Town□Dongguan City□Grina□China