

Trade Term: FOB 30 days/100%

Quantity: 1,000,000 Pcs

Payment: T/T 100% in advance

Lead Time: 30 days

Product Size: 37mm x 15mm x 20mm

Certification: UL, CE, ROHS, FCC, UKCA

Shipping: Fedex, DHL, UPS

MOQ: 500 PCS

Notes:

- Product Name: TENS 7000, Nicwell Care, TENS 3000
- Material: ABS
- TENS Electrode: DC 2.35mm, 2.0mm, 1150mm
- Brand: TENS, EMPI, Prospera
- Application: Pain Relief

Product Name: TENS Electrode



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



ISO13485

FCC

CE

ROHS



UKCA

EPR

EPR

ISO9001

--	--	--	--	--	--

Q1:

A: Heymed

Q2:

A: PayPal T/T L/C

Q3:

A: □□□□□□□□□□□□□□QC□□□□□□□□□□□□□□□□□□Rohs□SGS□CE
FCC□ISO9001:2008□□□□□□□□□□

Q4

A: 100 OEM/ODM

Q5: □□□□□?□□□□□□□□□□□□□□□□□□□□?

A: 2 4 12 DHL Fedex UPS TNT EMS

Q6: □□□□□□□□□□? □□□□□□□□□□□□□□□□□□?

A: