

8 Pin Lemo

Pin	8 Pin Lemo Connector - 6Pin Connector
Pin	8 Pin Lemo Connector 6 Pin Connector
Length	30cm
Material	TPU
Application	EEG ECG EMG
MOQ	100

8 Pin Lemo Connector 6 Pin Connector [8 Pin Lemo Connector](#) ECG EEG EKG EMG Lemo Connector 8 Pin Connector TPU TPE

DIN Connector [Lemo Connector - 40mm](#)

:





□□□□□□



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



