

品牌	10.0mm ISO13485 ISO9001 Dongguan China Factory
规格	
尺寸	10cm 20cm 50cm 100cm
材料	
应用	
包装	

10.0mm EMS/TENS 10.0mm

ISO13485 ISO9001 guan TPU PVC

LEMO 3.5mm 4.0 ECG

B2B OEM/ODM

•

- 10.0mm
-
-
- ISO13485 | ISO9001
- OEM/ODM MOQ





About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10years of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



ISO 13485

FCC

CE

ROHS



UKCA

EPR

EPR

ISO 9001

