

Trade Term: FOB 30 days

Quantity: 1,000,000 pcs

Payment: T/T

Lead Time: 30 days

Delivery Time: 3-7 days after 15-20 days

Certification: UL, CE, ROHS, FCC, UKCA

Logistics: Fedex, DHL, UPS

MOQ: 500 PCS

3.5mm 4.0mm EMS ECG **3.5 4.0mm**

2.0

12 RoHS Reach

ISO 13485 ISO 9001 ISO 14001 UL CQC CE SGS



FCC ISO9001:2008

Q4:

A: 100 OEM/ODM

Q5:

A: 2 4 12 DHL Fedex UPS TNT EMS

Q6:

A: