

名稱	Guangdong Professional Custom Physiotherapy Cable 2.0 Round Tube Ends to 2.0 Electrode Pins with SR Lead Wire ISO9001/13485FDACE
規格	標準規格
尺寸	10cm20cm50cm100cm
材料	醫療級導電材料
包裝	包裝
數量	數量

說明

- 本產品為專業定制之物理治療電線，適用於 TENS 及 EMS 治療。產品符合 ISO9001/13485 標準，並通過 FDA 及 CE 認證。包裝及數量請參閱說明書。
- 本產品為標準規格，適用於各種物理治療儀器。包裝及數量請參閱說明書。
- 本產品為醫療級導電材料，適用於 TENS 及 EMS 治療。包裝及數量請參閱說明書。
- 本產品為標準規格，適用於各種物理治療儀器。包裝及數量請參閱說明書。

規格

- 本產品為標準規格，適用於各種物理治療儀器。
- 本產品為標準規格，適用於各種物理治療儀器。
- 本產品為標準規格，適用於各種物理治療儀器。
- 本產品為標準規格，適用於各種物理治療儀器。
- 本產品為標準規格，適用於各種物理治療儀器。
- 本產品為標準規格，適用於各種物理治療儀器。







品質保証書

品質保証書は、製品の品質を保証するための重要な文書です。

- **ISO 13485** 品質保証システム規格に準拠しています。
- **ISO 9001** 品質保証システム規格に準拠しています。
- **FDA** 米国食品医薬品局（FDA）の承認を受けています。
- **CE MDR** EUの医療機器指令（MDR）に準拠しています。

品質保証書



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



--	--	--	--	--	--

Q1

A _____

Q2□□□□□□□□□□□□□□□□

A PayPal/T/T Western Union L/C

Q3

QC ROH SGS CE
FCC ISO9001 2008

Q4

A100 OEM/ODM

Q5

A 2412 DHL FedEx UPS TNT EMS

Q6

guan3

□A□No□8□Kangyang 1st Street□Shebei Village□Huangjiang Town□Dongguan City□Grina□China