

品名	DB15 13pin 4.0mm 中国东莞 EMS TENS
规格	
尺寸	10cm 20cm 50cm 100cm
颜色	
材质	
包装	
备注	

DB15 13pin 4.0mm SNAP 中国东莞 DB15 13pin 4.0mm 中国东莞 13pin 4.0mm 中国东莞

TPU TENS EMS ECG/EKG 2 ECG/EKG USB-C 4.0mm C

Dongguan OEM/ODM ISO13485 ISO9001



.



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



Q1
A

Q2
A PayPal T/T Western Union L/C

Q3
 QC ROH SGS CE
FCC ISO9001 2008

Q4
A 100 OEM/ODM

Q5
A 2 4 12 DHL FedEx UPS TNT EMS

Q6
 guan 3
 A No 8 Kangyang 1st Street Shebei Village Huangjiang Town Dongguan City Grina China