


規格	M12 M16 4 PIN 4.0mm ISO13485 Dongguan
材質	TPU / PVC
長度	10 cm 20 cm 50 cm 100 cm
顏色	白色 / 黑色
應用	醫療 / 工業
認證	ISO13485 / ROHS
包裝	1000PCS / 500PCS

Dongguan 4 M12 M16 4.0mm 2.0 3.5mm/4.0mm TENS EMS

S TPU PVC M12/M16 4 - Pin TEN EMS ECG EMG LEMO 4.0mm ECG

OEM ODM ISO13485 ROHS





About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



ISO13485

FCC

CE

RoHS



UKCA

EPR

EPR

ISO9001

