


Ag/AgCl Electrode Array for ECG, EEG, and EMG Applications

High-quality, disposable Ag/AgCl electrodes for medical monitoring. Available in various sizes and configurations. [View Details](#) [Contact Us](#)

Manufactured by Dongguan Factory, B2B OEM/ODM services available.






1. 產品名稱

2. 產品用途

- **ISO 13485** 產品符合 ISO 13485 標準
- **ISO 9001** 產品符合 ISO 9001 標準
- **FDA** 產品符合 FDA 標準
- **CE MDR** 產品符合 EU MDR 標準

3. 產品規格



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



ISO13485

FCC

CE

RoHS



UKCA

EPR

EPR

ISO9001

Q1
A

Q2
A PayPal T/T Western Union L/C

Q3
 QC ROH SGS CE
FCC ISO9001 2008

Q4
A 100 OEM/ODM

Q5
A 2 4 12 DHL FedEx UPS TNT EMS

Q6
 guan 3
 A No 8 Kangyang 1st Street Shebei Village Huangjiang Town Dongguan City Grina China