

Trade Term: FOB 30/60/90

Quantity: 1,000,000 Pcs

Payment: T/T

Lead Time: 30 Days

Production Capacity: 300,000 Pcs/Month

Certifications: UL, CE, ROHS, FCC, UKCA

Logistics: Fedex, DHL, UPS

MOQ: 500 PCS

Notes:

1. The price is based on 90 days payment term. If you want to pay 30 days, please contact us for the discount.
2. The price is based on the standard configuration. If you want to customize, please contact us for the quotation.
3. The price is based on the standard configuration. If you want to customize, please contact us for the quotation. 10% discount for 10% order quantity.

Remarks:


[90 - DC2.35 - 2.0 ECG 2.0 4](#)

12 Months Warranty, RoHS, Reach

ISO 13485, ISO 9001, ISO 14001, UL, CQC, CE, SGS



型号:DC3.5



--	--	--	--	--	--

Q1:

A: Heymed

Q2:

A: PayPal T/T L/C

Q3:

A: □□□□□□□□□□□□□□□□QC□□□□□□□□□□□□□□□□□□□□□□Rohs□SGS□CE

FCC ISO9001:2008

Q4:

A: 100 OEM/ODM

Q5:

A: 2 4 12 DHL Fedex UPS TNT EMS

Q6:

A: