

Trade Term: FOB

Quantity: 1,000,000

Payment: T/T

Lead Time: 30 days

Product Dimensions: 37x15x20mm

Certifications: UL, CE, ROHS, FCC, UKCA

Shipping: Fedex, DHL, UPS

MOQ: 500PCS

Product

- Heymed 2.0mm USB - TENS TENS 7000, Nicwell Care, TENS 3000
- Applications: TENS, EMS, E-stim, etc.
- TENS Applications: [TENS](#), TENS, TENS
- Tens Applications: TENS, EMS, E-stim, etc.
- Applications: [Tens](#), TENS, EMS, E-stim, etc.

Product [TENS](#) Applications:



Product

[3 Applications](#)

[USB C - 2.0 mm Applications](#)

Heymed Technology Applications [Applications](#) 3.5mm Applications Applications Heymed Applications



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



ISO13485

FCC

CE

ROHS



UKCA

EPR

EPR

ISO9001

