

產品規格:

產品名稱	DC5521 2 in 1 經皮神經電刺激器 - DC5521 2 in 1 經皮神經電刺激器
型號	DC5521
尺寸	30cm
重量	約 100g
功能	TENS EMS 經皮神經電刺激 ECG/EKG/EEG/EMG 監測 記錄
電池壽命	500h

特點:

- 產品: DC5521 2 in 1 經皮神經電刺激器 - DC5521 2 in 1 經皮神經電刺激器
- 12 個通道: 12 個通道, 每個通道可獨立設定刺激參數, 包括刺激強度、頻率、脈寬、刺激時間等。
- 多種刺激模式: 提供多種刺激模式, 包括 TENS、EMS、ECG、EKG、EEG、EMG 等。
- 數據記錄: 可記錄刺激數據, 包括刺激強度、頻率、脈寬、刺激時間等, 並可导出數據文件。
- 電池壽命: 電池壽命可達 500 小時。
- 尺寸: 尺寸為 30cm。
- 重量: 重量約 100g。
- 兼容性: 兼容 DC5521 2 in 1 經皮神經電刺激器。


規格:

- 產品名稱: DC5521 2 in 1 經皮神經電刺激器 - DC5521 2 in 1 經皮神經電刺激器
- 型號: DC5521
- 尺寸: 30cm
- 重量: 約 100g

認證:

- CE 認證: CE 認證 ECG 監測器
- FDA 認證: FDA 認證 (FDA 註冊) 經皮神經電刺激器
- RoHS 認證: RoHS 認證 (RoHS) 經皮神經電刺激器


包裝:







□□□□□□



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



