





00 00. 00



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



FAQs

Q1: What is Heymed?

A: Heymed is a professional medical device manufacturer and exporter. We have been in the industry for over 10 years. Our products are widely used in hospitals and clinics.

Q2: What is your payment term?

A: T/T, L/C, Western Union, etc.

Q3: What is your quality control system?

A: We have a strict quality control system. All our products are tested before shipment. We are ISO9001:2008 certified. We also have RoHS, SGS, CE, FCC, etc. certificates.

Q4: What is your lead time?

A: 100% on time. We are OEM/ODM manufacturer. We can produce 10000 units per month.

Q5: What is your shipping method?

A: We can ship by DHL, Fedex, UPS, TNT, EMS, etc. It takes 2-4 days to reach your door.

Q6: What is your price?

A: We are a professional manufacturer. Our price is very competitive. We can provide 3rd party inspection. We can also provide 100% AQL.