

ECG EKG EEG EMG

Model	ECG EKG EEG EMG - 2.54mm Dupont 1x1p 30cm Cable
Pin	2.54mm 1x1p 3.5mm/4.0mm 1x1p ecg cable
Cable Length	30cm 1x1p 30cm
Material	1x1p 30cm
Application	TENS EMS 1x1p/1x1p ECG/EKG/EEG/EMG 1x1p 30cm
MOQ	100PCS

ISO13485 Medical Device Manufacturer. [ECG EKG EEG EMG 1x1p 30cm](#) .


2.54mm 1x1p 30cm 1x1p 30cm 4.0mm ecg cable 1x1p 30cm 1x1p 30cm, 1x1p 30cm 4.0mm 1x1p 30cm, 1x1p 30cm d2.54mm 1x1p 30cm. 1x1p, 1x1p OEM ODM 1x1p 1x1p 1x1p 1x1p 1x1p [RJ45 - 1x1p 1x1p](#). 1x1p, 1x1p 1x1p 1x1p 1x1p 1x1p 1x1p. [1x1p 1x1p ECG 1x1p 1x1p 1x1p M5 1x1p 1x1p](#) 1x1p 1x1p 1x1p 1x1p 1x1p 1x1p 1x1p 1x1p 1x1p 1x1p.

1x1p 1x1p:





□□ □□. □□



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



歡迎來電 查詢

Q1: 貴公司 有 無 專 門 的 售 後 服 務 嗎 ？

A: 當然有。Heymed 有 專 門 的 售 後 服 務 團 隊， 專 門 為 各 個 客 戶 提 供 售 後 服 務。 如 有 任 何 問 題， 請 即 與 我 們 聯 繫。

Q2: 貴 公 司 有 無 售 後 服 務 ？

A: 當然， T/T, 信 用 證， 或 者 現 款 交 貨 L/C

Q3: 貴 公 司 有 無 售 後 服 務 團 隊 ？

A: 當然 有。 我 們 有 專 門 的 售 後 服 務 團 隊， 專 門 為 各 個 客 戶 提 供 售 後 服 務。 我 們 通 過 Rohs, SGS, CE FCC, ISO9001:2008 等 證 書 來 證 明 我 們 的 質 量。

Q4: 貴 公 司 有 無 售 後 服 務 ？

A: 當然有， 我 們 100% 專 門 為 各 個 客 戶 提 供 售 後 服 務。 我 們 有 OEM/ODM 等 售 後 服 務。 如 有 任 何 問 題， 請 即 與 我 們 聯 繫。

Q5: 貴 公 司 有 無 售 後 服 務 ？

A: 我 們 有 專 門 的 售 後 服 務 團 隊， 專 門 為 各 個 客 戶 提 供 售 後 服 務。 我 們 通 過 DHL, Fedex, UPS, TNT, EMS 等 證 書 來 證 明 我 們 的 質 量。

Q6: 貴 公 司 有 無 售 後 服 務 ？

A: 當然有。 我 們 有 專 門 的 售 後 服 務 團 隊， 專 門 為 各 個 客 戶 提 供 售 後 服 務。 我 們 通 過 3rd 等 證 書 來 證 明 我 們 的 質 量。 我 們 有 專 門 的 售 後 服 務 團 隊， 專 門 為 各 個 客 戶 提 供 售 後 服 務。 我 們 通 過 100% 專 門 為 各 個 客 戶 提 供 售 後 服 務。