

ECG 3.5mm/4.0mm/5.0mm

ECG 3.5mm/4.0mm/5.0mm	ECG 3.5mm/4.0mm/5.0mm
30cm	30cm
TENS EMS	TENS EMS
MOQ	100PCS

ECG 3.5mm/4.0mm/5.0mm 30cm TENS EMS 100PCS. OEM ODM EEG ECG. MOQ: 100PCS.

A: 4.0

B: ECG

MOQ:





□□ □□. □□



## About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M<sup>2</sup>  
Factory space with 10years of production experience



100+  
Employees



OEM/ODM  
Support cable customized



TOP 500  
Bands approved manufacturer



## Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



ISO 13485

FCC

CE

ROHS



UKCA

EPR

EPR

ISO 9001

□□□□ □□

**Q1:** □ □□□□ □□ □□ □□□ □ □□□□?

A: □□□□. Heymed□ □□□ □□ □□ □□ □□□□□, □□□□ □□ □□ □□□□ □□□. □□□□ □□ □□□ □□□ □□□□□.

**Q2:** □□ □□ □□□ □□□□□?

A: □□□, T/T, □□ □□, □□□ □□□ □ L/C

**Q3:** □□□ □□ □□□ □□□□ □□□ □□ □□□□?

A: □□□ □□□□□. □□□□ □□ □□ □□□□ □□ QC □□□□ □□□□ □□ □□□ □□ □□□□ □□□□□. □□ □□□ Rohs, SGS, CE FCC, ISO9001:2008 □□□ □□□□□□. □

**Q4:** □□? □! □□□ □□□ □□□□□ □□□□ □□ □□ □□□ □□□ □□□ □□ □□□□?

A: □□□□□, □□□ 100□ □□ □□□ □□ □□□□□. □□□ OEM/ODM □□□□ □□□□□. □□□ □□□ □□ □□ □□□□□ □□□□□□□.

**Q5:** □□ □□ □ □□? □ □□□ □□ □□ □□□□ □□□□□?

A: □□ □□□ □ 2-4□□□□ □□ □ 12□□ □□□ □□□ □ □□□□. □□□ DHL, Fedex, UPS, TNT, EMS□ □□□□ □□□□□□.

**Q6:** □□□ □□□ □□□□? □□□ □□□ □□□ □ □□□□?

A: □□□□□. □□□ □□ □□ □□□□ □□ □□□ □□□ □□ □□□□. □□ □□□ 3rd□ □□□□□. □□ □□□ □□□ □□□ □□□□ □□1□□ 8□ A□ 1□