

[illegible]

EEG 是研究腦電活動的一種方法，它記錄了大腦皮層神經元的電活動。EEG 的優點是時間分辨率高，可以捕捉到毫秒級的電活動變化。然而，EEG 的空間分辨率較低，難以精確定位腦內活動的具體位置。因此，EEG 通常與其他成像技術（如 fMRI、PET）結合使用，以實現更全面的腦功能研究。

本公司為 100% 中國公司，所有產品均由中國製造。Dongguan 工廠位於中國東莞，所有產品均由中國製造。本公司為 OEM/ODM 工廠。





□□ □ □□ □□ :

□ □□□□ □□ □□□□ □□ □□ □ □□ □□□□ □□ □□□ □□□□□□□□.

- **ISO 13485** : □□□ □□ □□ □□ □□ □□□□ □□ □□ □□□□□□.
- **ISO 9001** : □□□ □□ □□ □□□□□ □□□□□□□□□□.
- **FDA** □□ : □ □□□ □□ □□□ □□□ □□□□ □□ □□ □□ □□□□□□.
- **CE MDR** □□ : □□ □□□ □□□□ □□□ □□□□ □□ EU □□ □□ □□ (MDR)□ □□□□□□.

□□ □□. □□



About Us

10+ Years Experienced Professional Manufacturer
of Customized Cables.



2000 M²
Factory space with 10 years
of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved
manufacturer



Certificate

We followed international ISO management Standard to make our custom cables,
medical cables and chargers with good quality and best price.



FAQ

Q1 : □ □□□ □□ □□ □□□ □□□□?

A : □□, Heymed□ □□□ □□□ □□ □□ □□ □□□□. □□□ □□□□ □□ □□□ □□□□□□. □□□ □□□ □□□ □□□ □□□□□□.

Q2 : □□ □□ □□□ □□□□□□□?

A : PayPal, T/T, ☐ ☐, Western Union ☐ L/C.

Q3 : □□□ □□ □□ □□□□ □□□□□□□□?

□ : □□□ □□□□□. □□ □□ □□ □□□□□□ QC □□□ □□□□ □□□ □□ □□□ □□□□□. □□□ ROH, SGS, CE FCC, ISO9001 : 2008 □□□ □□□□□. □

[illegible]

A : 00000. 000 000 100 0 000 00000. OEM/ODM 0000 00000. 000 000 00 00 00 00000 000000.

Q5 : ພວກເຮົາໄດ້ສ້າງຕັ້ງຄອບຄົວຂຶ້ນມາແລ້ວ ແຕ່ພວກເຮົາຍັງບໍ່ທັນມີຄວາມສະຫງົບໃຈ ເປັນຫຍັງ?

A : 2-4 12 DHL, FedEx, UPS, TNT, EMS

Q6 : ?

□ : □□. □□□ Dongguan□□ □□□ □□ □□□ □□ □□□□. □□ □□□ □ □□□ □□□□. A, No. 8, Kangyang 1st Street, Shebyi Village, Huangjiang Town, Dongguan City, Guangdong Province, □□