





□□ □□. □□



About Us

10+ Years Experienced Professional Manufacturer of Customized Cables.



2000 M²
Factory space with 10years of production experience



100+
Employees



OEM/ODM
Support cable customized



TOP 500
Bands approved manufacturer



Certificate

We followed international ISO management Standard to make our custom cables, medical cables and chargers with good quality and best price.



ISO 13485

FCC

CE

ROHS



UKCA

EPR

EPR

ISO 9001

FAQ

Q1: 貴公司はどのような製品を扱っていますか？

A: はい、Heymedの製品は、主に医療機器、美容機器、産業機器の分野で活躍しています。品質が高く、信頼性の高い製品を提供しています。

Q2: 貴社の支払い方法はありますか？

A: PayPal, T/T, 銀行振込, 信用状 (L/C) などの支払い方法に対応しています。

Q3: 貴社の製品はどのような規格に準拠していますか？

はい、製品は厳格な品質管理システム (QC) を実施しています。また、RoHS, SGS, CE, FCC, ISO9001:2008 などの国際規格に準拠しています。

Q4: 貴社の MOQ (Minimum Order Quantity) はどれくらいですか？

A: 通常、100個以上の MOQ を設定しています。OEM/ODM 対応も可能です。お問い合わせください。

Q5: 貴社の製品のリードタイムはどれくらいですか？

A: 通常、2~4週間、最大12週間以内のリードタイムです。DHL, Fedex, UPS, TNT, EMS などの物流サービスも提供しています。

Q6: 貴社の製品はどのような保証期間がありますか？

はい、製品は1年間の保証期間を提供しています。また、アフターサービスも充実しています。お問い合わせください。